

STOMADENT TU CART

Mobile dentist's table with bottom-led hoses equipped with up to 5 instruments (3 of them optic)




Pedal
foot
control

Basic version includes:

- 3-way syringe
- Midwest hose for optic turbine with contacts
- Foil keyboard control (Model S200 only)
- Built-in source of pressure water for instrument cooling (bottle)
- Pedal foot control

Additional optional instruments and functions can be added, according to customer-requests.



Model
S100



Technical data:

Power supply: 230 V \pm 10%
Frequency: 50 Hz \pm 2%
Maximum power input at 230V/50Hz: 50 VA + 10%
Unit weight (with/without compressor): 49 kg/16 kg
Type of shock protection: Class I Equipment
Degree of shock protection: Applied parts of B type
Recommended max. load of table: 1,5 kg
Input air pressure: 0,45 - 0,8 MPa

STOMADENT TU CART is available with or without integrated compressor.

Model
S200

